

# gp120(MN)



---

Product Name	gp120(MN)(HIV-1/Clade B)
Catalog No	IT-001-002MNp
Source	Purified protein from transfected 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 Clade B gp120(MN) protein (a.a.30-513 ) (Genebank accession no. AF075722)
Applications	Antibody ELISA, WB standard, , immunogen, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

