




HA(Δ TM) (H₁₅N₈)

Product Name	HA(Δ TM) (H ₁₅ N ₈)
Catalog No	IT-003-016 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM) (H ₁₅ N ₈) (A/duck/AUS/341/1983(H15N8)) protein (a.a.18-534) (Genebank No.ABB88132)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM) (H₁₅N₈) protein

Note: