



HA(Δ TM) (H₁₆N₃)

Product Name	HA(Δ TM) (H ₁₆ N ₃)
Catalog No	IT-003-017 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM) (H ₁₆ N ₃) (A/black-headed gull/Sweden/5/99(H16N3)) protein (a.a.16-530) (Genebank No.AAV91217)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

