



# HA( $\Delta$ TM) (H<sub>3</sub>N<sub>2</sub>)

---

Product Name	HA( $\Delta$ TM) (H <sub>3</sub> N <sub>2</sub> ) (A/Brisbane/10/2007)
Catalog No	IT-003-0042 $\Delta$ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM) (H <sub>3</sub> N <sub>2</sub> ) (A/Brisbane/10/2007) protein (a.a.17-530) (Genebank No. EU199248)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM) (H<sub>3</sub>N<sub>2</sub>) (A/Brisbane/10/2007) protein

Note: