




HA(Δ TM)(H₅N₈)

Product Name	HA(Δ TM)(H ₅ N ₈)
Catalog No	IT-003-020 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(H ₅ N ₈) (A/duck/NY/191255-59/02(H5N8)) protein (a.a.1-490) (Genebank No. AAP72011)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₅N₈) protein

Note: