




HA(Δ TM)(H₅N₂)

Product Name	HA(Δ TM)(H ₅ N ₂)
Catalog No	IT-003-018 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(H ₅ N ₂) (A/American green-winged teal/California/HKWF609/2007(H5N2)) protein (a.a.19-506) (Genebank No. ACF47563)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₅N₂)

Note: