



HA(Δ TM)(H₅N₃)

Product Name	HA(Δ TM)(H ₅ N ₃)
Catalog No	IT-003-019 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(H ₅ N ₃) (ectodomain of HA) (A/duck/Hokkaido/167/2007(H5N3)) protein (a.a.18-500) (Genebank No. BAG07130)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

