



# HA( $\Delta$ TM)(A/2007(H<sub>6</sub>N<sub>1</sub>))

---

Product Name	HA( $\Delta$ TM )(A/northern shoveler/California/HKWF115/2007(H <sub>6</sub> N <sub>1</sub> ))
Catalog No	IT-003-021 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM )(A/northern shoveler/California/HKWF115/2007(H <sub>6</sub> N <sub>1</sub> )) protein (a.a. 17-529) (Genebank No. ACE81692).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)( A/2007(H<sub>6</sub>N<sub>1</sub>))