

HA(Δ TM)(A/Hongkong/1073/1999(H₉N₂))



Product Name	HA(Δ TM)(A/Hongkong/1073/1999(H ₉ N ₂))
Catalog No	IT-003-0091 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged H ₉ (Δ TM)(A/Hongkong/1073/1999(H ₉ N ₂)) protein (a.a. 19-530)(Genebank No. NP_859037).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

