

HA(Δ TM)(A/Chicken/Hong Kong/G9/97(H₉N₂))



Product Name	HA(Δ TM)(A/Chicken/Hong Kong/G9/97(H ₉ N ₂))
Catalog No	IT-003-0092 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(Δ TM)(A/Chicken/Hong Kong/G9/97(H ₉ N ₂)) protein (a.a. 19-530) (Genebank No. AAF00701).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/Chicken/Hong Kong/G9/97(H₉N₂))