

H₁₀(ΔTM) (H₁₀N₇)



Product Name	H ₁₀ (ΔTM)(A/blue-winged teal/Louisiana/Sg-00073/07(H10N7))
Catalog No	IT-003-011ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged H ₁₀ (ΔTM) (A/blue-winged teal/Louisiana/Sg-00073/07(H10N7)) protein (a.a. 18-526) (Genebank No. ADH96304)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

