



H₁₁(ΔTM) (H₁₁N₂)

Product Name	H ₁₁ (ΔTM)(A/duck/Yangzhou/906/2002(H11N2))
Catalog No	IT-003-012ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged H ₁₁ (ΔTM)(A/duck/Yangzhou/906/2002 (H11N2)) protein (a.a. 17-528) (Genebank No. AAU85533)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E:H₁₁(ΔTM)(A/duck/Yangzhou/906/2002(H11N2))