

H₁₂(ΔTM) (H₁₂N₅)



Product Name	H ₁₂ (ΔTM)(A/green-winged teal/ALB/199/1991(H12N5))
Catalog No	IT-003-013ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged H ₁₂ (ΔTM)(A/green-winged teal/ALB /199/1991(H12N5)) protein (a.a. 17-530) (Genebank No. ABB88110)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

