

# HA<sub>1</sub>(H13N8)(A/black-headed gull/Netherlands/1/00)



Product Name	HA <sub>1</sub> (H13N8)(A/black-headed gull/Netherlands/1/00)
Catalog No	IT-003-014p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (H13N8)(A/black-headed gull/Netherlands/1/00) protein (a.a.17-341) (Genebank #: <a href="#">AAV91212</a> )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(H13N8)(A/black-headed gull/Netherlands/1/00)