

HA₁(H15N8)(A/duck/AUS/341/1983)



| | |
|----------------|---|
| Product Name | HA ₁ (H15N8)(A/duck/AUS/341/1983) |
| Catalog No | IT-003-016p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | C-terminal 6xHis tagged HA ₁ (H15N8)(A/duck/AUS/341/1983) protein (a.a.18-349) (Genebank #: ABB88132) |
| Applications | WB standard, antibody ELISA, immunogen, etc |
| Purity | >95% |
| Storage | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: HA₁(H15N8)(A/duck/AUS/341/1983)