

# H3( $\Delta$ TM)(H3N2)(A/Wuhan/359/95)



---

|                |  |
|----------------|--|
| Product Name   | H3( $\Delta$ TM)(H3N2)(A/Wuhan/359/95)   |
| Catalog No     | IT-003-0049 $\Delta$ TMp   |
| Source         | Recombinant protein purified from 293 cell culture   |
| Formulation    | Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS  |
| Protein        | C-terminal 6xHis tagged H3( $\Delta$ TM)(H3N2)(A/Wuhan/359/95) protein (a.a.17-529) (Genebank #: AAT64865)           |
| Applications   | WB standard, antibody ELISA, immunogen, etc  |
| Purity         | >95%   |
| Storage        | Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

