

HA(Δ TM) (H₉N₂) (A/Hong Kong/33982/2009)



Product Name	HA(Δ TM) (H ₉ N ₂)(A/Hong Kong/33982/2009)
Catalog No	IT-003-0094 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of HA(H ₉ N ₂)(A/Hong Kong/33982/2009) protein (a.a. 17-523) (Genebank #: EPI241740)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H₉N₂)(A/Hong Kong/33982/2009)