

HA(Δ TM)(A/NewJersey/8/1976)(H1N1)



Product Name	HA(Δ TM)(A/NewJersey/8/1976)(H1N1)
Catalog No	IT-003-SW5 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/NewJersey/8/1976)(H1N1) protein (a.a.18-529) (Genebank No. ACQ99821).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(A/NewJersey/8/1976)(H1N1)