




S (Beijing02)(SARS-CoV)

Product Name	S (1-1190)(Beijing02)(SARS-CoV)
Catalog No	IT-002-001p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged Spike protein (amino acid 1-1190) (SARS-CoV) (Beijing02) (Genebank No.AY278487)
Purity	≥ 95%
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified S (1-1190)(SARS-CoV) protein.