



gp120(CAP45.2.00)(HIV-1)

Product Name	gp120(CAP45.2.00)
Catalog No	IT-001-RC11p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 (CAP45.2.00)(Clade C) protein (a.a.34-518) (Genebank No. DQ435682).
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(CAP45.2.00) protein

Note: