




# gp120 (HIV-1/Clade C)

---

Product Name	gp120 (HIV-1/Clade C)
Catalog No	IT-001-003p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 0.1% BSA and 25% glycerol.
Protein	6xHis tagged HIV-1 gp120 protein (a.a. 24-515) of strain C.ZA.1197MB from South Africa (Genebank #AY463234)
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified gp120 (Clade C) protein.