

# Rabbit IgG



---

Product Name	Rabbit IgG
Catalog No	IT-000-003
Source	SPF Rabbit(New Zealand White)
Isotype	IgG
Formulation	Each vial contains 1 mg purified IgG in 0.5 ml PBS
Applications	Negative control for ELISA, Neut., IP, IF etc.
Storage	Store at 4°C; stable for 3 months from the date of shipment. For long term storage, aliquot and freeze at -80°C. <b>Avoid repeated freeze and thaw.</b> Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.