



# Human RBP4

---

Product Name	Human retinol binding protein 4
Catalog No	IT-010-001p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged hRBP4 protein (a.a.19-201) (Genebank No. NM_006744)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

