



Human BNP

Product Name	hBNP (preproprotein)
Catalog No	IT-010-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged human natriuretic peptide precursor B proprotein protein (a.a.27-134) (Genebank No. NP_002512).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: human BNP protein (a.a.27-134)

Note: