



Anti- Human RBP4

Product Name	Anti- Human RBP4
Catalog No	IT-010-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Human retinol binding protein 4(hRBP4) (Genebank accession # NM_006744)
Applications	ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with Human retinol binding protein 4(hRBP4). Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.