



# Anti- hBNP (preproprotein)

---

Product Name	Anti- Human BNP (preproprotein)
Catalog No	IT-010-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Human natriuretic peptide precursor B proprotein protein (a.a.27-134) (Genebank No. <a href="#">NP_002512</a> ).
Applications	ELISA(1:1000-1:100,000), etc.
Specificity	Reacts with Human BNP (preproprotein) and mature products
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-010-002**

