

Human CD137 Ectodomain



Product Name	Human CD137 Ectodomain
Catalog No	IT-300-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged human CD137 ectodomain protein (a.a. 24-190), also called "4-1BB", which is the member of TNF receptor superfamily (Genebank Accession No. BC006196).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

