

Human CD3 Complex($\epsilon\gamma\delta$)



Product Name	Human CD3 Complex($\epsilon\gamma\delta$)
Catalog No	IT-300-hCD3p-egd
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of human CD3 Complex($\epsilon\gamma\delta$) proteins: ϵ (a.a. 20-126) (Genebank #: X03884) in complex with γ (a.a.28-113)(Genebank #: NP_000064), and δ (a.a. 24-103) (Genebank #:NP_000723).
Applications	WB standard, antibody ELISA, immunogen, T cell functional analysis, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human CD3 Complex($\epsilon\gamma\delta$)