

LecCapture™



Product Name	LecCapture™
Catalog No	IT-200-001
Source	Purified by affinity chromatography
Formulation	10 µg/ml Lectin in phosphate buffer, 100 ml in total.
Applications	Coating reagent for capturing glycoprotein for ELISA analysis. Use 100 µl/well for 96-well-dish.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Wear glove when handling.
Use Limitation	For research use only, not for use in diagnostic procedures.