




# Anti-HBsAg(Pre-S<sub>2</sub> Specific)

---

Product Name	Anti-HBsAg (Pre-S <sub>2</sub> Specific) (HBV)
Catalog No	IT-004-0011
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HBsAg Pre-S <sub>2</sub> protein (Genebank accession #AY040627)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with Pre-S <sub>2</sub> of HBsAg antigen from HBV. Not cross-react to HBsAg antigen.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;  
1: 293 cell expressing HBsAg  
2: 293 cell expressing Pre-S<sub>2</sub>

Note: three different forms of Pre-S<sub>2</sub> detected.