

Anti-NS5b (HCV)



Product Name	Anti- NS5b (HCV)
Catalog No	IT-004-022
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide from a conserved region of NS5b polypeptide of Hepatitis C virus (a.a. 2444-2454) (CSMSYSWTGAL).
Applications	WB(1:1000), ELISA(1:2000-1:10000), IP, etc.
Specificity	Reacts with NS5b polypeptide of Hepatitis C virus (a.a. 2444-2454) (CSMSYSWTGAL).
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.

ELISA Profile of IT-004-022

