

# Anti-Human TIM-3(Ectodomain)



Product Name	Anti-Human TIM-3 (Ectodomain)
Catalog No	IT-004-TIM-3M1
Source	Mouse monoclonal antibody, clone 2D9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing ectodomain of human TIM-3 protein (a.a. 22-200) (Genebank #: AAL65158).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with human TIM-3 protein. TIM-3 is the member of T-cell immunoglobulin and mucin domain-containing proteins. TIM-3 is also a putative receptor of HAV. Cross-reactivity not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-004-TIM-3M1

