




Anti-VP4 (HAV)

Product Name	Anti-VP4 (HAV)
Catalog No	IT-004-013
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hepatitis A virus VP4 protein (Genebank accession #AAM08224)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HAV VP4 and VP4 precursor proteins. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing VP4 (HAV)

Note: