




Anti-gG (HSV-1)

Product Name	Anti-gG (HSV-1)
Catalog No	IT-005-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Glycoprotein G of human herpes virus 1 (HSV-1) (Genebank No. CAF24800)
Applications	Western Blotting (1:100-1:1000)
Specificity	Reacts with glycoprotein G of human herpes virus 1 (HSV-1). Cross-reactivity to gG from other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing gG

Note: