

Anti-VP1(Norovirus/ Hu/Shanghai/ SH5/2009/CHN)



| | |
|----------------|--|
| Product Name | Anti-VP1(Norovirus/Hu/Shanghai/SH5/2009/CHN) |
| Catalog No | IT-013-002 |
| Source | Rabbit Polyclonal Antibody |
| Formulation | Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide. |
| Immunogen | DNA vaccine expressing VP1 of Norovirus (Hu/Shanghai/SH5/2009/CHN) protein(a.a.1-540) (Genebank accession # ADF47128). |
| Applications | ELISA(1:1000-1:100,000), WB, IP, Neutralizing assay, etc. |
| Specificity | Reacts with human Norovirus VP1 protein. Cross-reactivity to other family member not tested. |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

