

# Anti-VP1(MPV-2)



---

Product Name	Anti-VP1(MPV-2)
Catalog No	IT-015-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing VP1 of Mouse parvovirus 2 protein (a.a.1-729) (Genebank accession # YP_656491).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with human APOA1 protein. Cross-reactivity to other family member not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

