


E2 glycoprotein (Mayaro virus)



| | |
|----------------|--|
| Product Name | E2 glycoprotein (Mayaro virus) |
| Catalog No | IT-022-025Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water. |
| Protein | C-terminal V5-6xHis tagged E2 glycoprotein (Mayaro virus) antigen (a.a.1-364) (Genbank Accession No. NP_740693). |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Purity | >90% , based on SDS PAGE |
| Storage | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: E2 glycoprotein (Mayaro virus)