

NS (Rift Valley Fever Virus)



Product Name	NS (Rift Valley Fever Virus)
Catalog No	IT-017-003Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged non-structural protein [Rift Valley Fever Virus] protein (a.a. 1-265) (Genebank Accession No. ABD38738).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

