

Glycoprotein G1 (SFTS Virus)



Product Name	Glycoprotein G1 (SFTS Virus)
Catalog No	IT-017-004p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged Phlebovirus glycoprotein G1 (SFTS virus HN6) (a.a.20-452) (Genebank No. ADZ04490).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Glycoprotein G1 (SFTS Virus)