




NP (Schmallenberg virus)

Product Name	NP (Schmallenberg virus)
Catalog No	IT-017-010p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged NP (Schmallenberg virus) protein (a.a. 1-233) (see amino acid sequence below)(Genbank Accession No. CCF55031)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ

mssqfifedvpqrnaatfnpevgvafigkygqqlnfgvarvfflnqkkakmvlhktapsvdltfggvkftvnnhfpq
yvsnpvpdnaitlhrmsgylarwiadtckasvklaeasaqivmplaevkgctwadgytmylgfapgaemfldafdfypl
viemhrvlkdnmdvnfmkkvlrqrygtmtaewmtqkiteikaafnsvgqlawaksgfspartflqqfginigs