

VP2 protein (Canine parvovirus)



Product Name	VP2 protein (Canine parvovirus)
Catalog No	IT-015-012Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged VP2 protein [Canine parvovirus] protein (a.a. 1-586) (Genebank Accession No. ACS50333).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: VP2 protein (Canine parvovirus)