

P3C (HAV)



Product Name	P3C (HAV)
Catalog No	IT-004-054Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 10mM CBB, pH9.6, 0.1% SDS and 50% glycerol. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.Coli derived 40 kDa recombinant protein contains the P3C immunodominant regions, amino acids 1643-1743.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.