




gp120(SF162)(HIV-1)

Product Name	gp120(SF162)
Catalog No	IT-001-0028p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120(SF162)(Clade B) protein (a.a.34-518) (Genebank No. M65024).
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Immune
Technology