



# Gag(HIV-1/Clade C)

---

Product Name	Gag(HIV-1/Clade C)
Catalog No	IT-001-009Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged Gag (HIV-1/Clade C) protein (Genebank accession No. AAB61122 )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at -80°C; <u>AVOID REPEATED FREEZE/THAW</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Gag (HIV-1/Clade C) protein

Note: