




p17 (HIV-1/Clade C)

Product Name	p17 (HIV-1/Clade C)
Catalog No	IT-001-011Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS with 8M urea. Reconstitute the protein with 100µl of Millipore water
Protein	C-terminal 6xHis tagged p17 (HIV-1/Clade C) matrix protein (Genebank accession no. AY463234)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Store at -80°C; <u>AVOID FREEZING/THAWING</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p17 (HIV-1/Clade C) protein

Note: