




p24 protein (HIV-1/HXBc2)

Product Name	p24 (HIV-1/HXBc2)
Catalog No	IT-001-017p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged p24 protein (Genebank accession #K03455)
Purity	≥ 95%
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified p24(HXBc2) protein