




p27 (SIV/mac239)

Product Name	p27 (SIV/mac239)
Catalog No	IT-001-018Ep
Source	Recombinant viral protein purified E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in PBS. Add 0.1 ml of ultrapure H ₂ O to reconstitute the protein after receiving. Spin briefly before reconstitution.
Protein	C-terminal 6xHis tagged SIV/mac239 p27 protein (Genbank No. M33262)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	> 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p27 (SIV/mac239)

Note: