

p24 Capsid Protein (HIV-1/Clade A) (Q23-CxC)



Product Name	p24 Capsid Protein (HIV-1/Clade A)(Q23-CxC)
Catalog No	IT-001-021p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged p24 Capsid Protein (HIV-1/Clade A)(Q23-CxC) protein (Genebank #: AAC97528).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p24 Capsid Protein (HIV-1/Clade A)
(Q23-CxC)