

# gp140(SF162) (HIV-1/Clade B)



---

Product Name	gp140(SF162) (HIV-1)
Catalog No	IT-001-143p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140 (SF162) (HIV-1) protein (a.a. 30-671) (Genebank No. AAT67508)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

