

# gp120(HIV-1/Clade CRF02\_AG/ 01CM\_0002BBY)



---

Product Name	gp120(HIV-1/Clade CRF02_AG/ 01CM_0002BBY)
Catalog No	IT-001-154p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120(HIV-1/Clade CRF02_AG/ 01CM_0002BBY) protein (a.a.30-506) (Genebank #: AY736843)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(HIV-1/Clade CRF02\_AG/  
01CM\_0002BBY)